



# User Guide Addendum Smart Tip (3x6.9mm Sharp)



	mm (in)		
	A	B	C
STT0513	4.7 (0.19)	5.7 (0.22)	7.8 (0.31)

Dimensions are based on 35-second cooling cycle performed in agarose gel.

This document is provided as a supplement to the User Guide. Always refer to the User Guide for complete operating instructions, warnings and precautions.



**WARNING!** To reduce likelihood of skin injury, skin warmer must be pressed against the skin when the needle is inserted before running a treatment cycle.

### Procedure:

- Prep cycle required for this tip type.